

Bone&JointScience

Our Innovation in Focus

Fatigue testing results of the newly redesigned PROMOS[◇] STANDARD Modular Shoulder System

Chris Moore¹, Jacob Cartner¹

¹ Smith & Nephew, Inc., Memphis, TN, USA

1 Research goal

To determine if a redesigned total shoulder arthroplasty (TSA) system can mitigate the risk of humeral head instability due to screw loosening and fracture.

2 Type of evidence



3 Clinical relevance

- The original PROMOS[◇] stem construct (Smith & Nephew, Inc.) has demonstrated acceptable clinical performance and return to function.¹ However, screw loosening and fracture are documented complications that can lead to humeral head instability and procedure failure in some instances.
- A newly developed design (Figure 1) was tested under worst case loading conditions to fully assess component fatigue performance and screw fracture risk.

4 Key result

- No fractures were observed for the new design following 7.5 million testing cycles, simulating 15-years of clinical use. Furthermore, there was no loss of screw torque.
- These results suggest that screw loosening and fracture are not expected to cause humeral head instability and failure in vivo.

5 Important considerations

- Additional research would be required to confirm component performance during mid to long-term clinical follow-up.



Figure 1: PROMOS[◇] STANDARD Modular Shoulder System (redesigned inclination set circled).

Background


The PROMOS° STANDARD Modular Shoulder System (Smith & Nephew, Inc., Memphis, TN, USA) is a shoulder arthroplasty implant featuring an adjustable inclination set and a clinically successful diaphyseal, rectangular cementless stem fixation design.² Biomechanical testing has demonstrated that this design has a higher torsional stability than the cylindrical shafts, without increased fracture risk.³ Clinically, the PROMOS stem provides intraoperative flexibility when addressing native anatomy, and has demonstrated good functional outcomes.¹ However, an analysis of complaint data and failure modes for the original stem construct revealed cases where loosening and slipping of the saddle joint caused failure of the inclination set-screw. This can lead to instability of the humeral head and clinical failure. Therefore, the PROMOS inclination set was redesigned for improved stability and fatigue strength (Figure 2). Construct changes include:

- Removal of transverse screw hole and polyethylene pin.
- Removal of saddle joint.
- Addition of double cone taper.
- Screw is now manufactured with rolled threads, as opposed to weaker, machined threads.

The purpose of the current study was to assess the mechanical strength of the re-designed system under simulated worst-case clinical loading conditions.

Methods

- A published in vitro model of the shoulder was referenced so that a rigorous and mechanically appropriate testing load could be applied to the newly redesigned components.⁴
- A total of six PROMOS STANDARD stem constructs were tested.
- Testing set-up generated a loosening torque about the longitudinal axis of the screw. Simulation protocol was designed to simulate 15-years of clinical use (maximum of 7.5 million joint cycles).
- All components were tested in a lubricated, physiologically correct liquid environment set at human body temperature (37°C).
- Set screw torque was assessed using a torquemeter before (5.6 N-m set point) and after testing to assess loosening. After the completion of testing, all modular components were disassociated and imaged to assess interface integrity.



Results

- No failures were observed during construct fatigue testing (Figure 3).
- The maximum testing cycles achieved was 7.5 million, designed to simulate 15-years of clinical use.
- Screw torque was largely unchanged following completion of post-fatigue testing ($5.54 \text{ N m} \pm 0.52 \text{ N m}$; $p > 0.05$).
- Taper disassociation was successfully achieved for each construct after completion of testing. Furthermore, imaging results of tapered surfaces did not show notable patterns of damage or wear.




Figure 3: Fatigue testing results for PROMOS° STANDARD stem constructs. All stem constructs (a–f) reached the maximum number of cycles (7.5 million, shown by green line) without failure.

Conclusion

The current testing results suggest that screw loosening and fracture are not expected to occur, and therefore are not expected to contribute to humeral head instability and clinical failure following implantation of the redesigned PROMOS° STANDARD Modular Shoulder System. However, additional research is necessary to confirm this finding over mid to long-term clinical follow-up.

References

1. Schumann K, Flury MP, Schwyzer HK, Simmen BR, Drerup S, Goldhahn J. Sports activity after anatomical total shoulder arthroplasty. *Am J Sports Med.* 2010;38(10):2097–2105.
2. Delaunay CP, Kapandji AI. Primary total hip arthroplasty with the Karl Zweymuller first-generation cementless prosthesis. A 5- to 9-year retrospective study. *J Arthroplasty.* 1996;11(6):643–652.
3. Flury MP, Schmoelz W, Schreiber U, Goldhahn J. Biomechanical testing of rectangular humeral shaft prosthesis: higher torsional stability without increased fracture risk. *Arch Orthop Trauma Surg.* 2011;131(2):267–273.
4. Bergmann G, Graichen F, Bender A, Kaab M, Rohmann A, Westerhoff P. In vivo glenohumeral contact forces--measurements in the first patient 7 months postoperatively. *J Biomechanics.* 2007;40(10):2139–2149.

Disclaimer Great care has been taken to maintain the accuracy of the information contained in the publication. However, neither Smith & Nephew, nor the authors can be held responsible for errors or any consequences arising from the use of the information contained in this publication. The statements or opinions contained in editorials and articles in this journal are solely those of the authors thereof and not of Smith & Nephew. The products, procedures, and therapies described are only to be applied by certified and trained medical professionals in environments specially designed for such procedures. No suggested test or procedure should be carried out unless, in the reader's professional judgment, its risk is justified. Because of rapid advances in the medical sciences, we recommend that independent verification of diagnosis, drugs dosages, and operating methods should be made before any action is taken. Although all advertising material is expected to conform to ethical (medical) standards, inclusion in this publication does not constitute a guarantee or endorsement of the quality or value of such product or of the claims made of it by its manufacturer. Some of the products, names, instruments, treatments, logos, designs, etc. referred to in this journal are also protected by patents and trademarks or by other intellectual property protection laws even though specific reference to this fact is not always made in the text. Therefore, the appearance of a name, instrument, etc. without designation as proprietary is not to be construed as a representation by the publisher that it is in the public domain. This publication, including all parts thereof, is legally protected by copyright. Any use, exploitation or commercialization outside the narrow limits of copyrights legislation, without the publisher's consent, is illegal and liable to prosecution. This applies in particular to photostat reproduction, copying, scanning or duplication of any kind, translating, preparation of microfilms and electronic data processing and storage. Institutions' subscriptions allow to reproduce tables of content or prepare lists of articles including abstracts for internal circulation within the institutions concerned. Permission of the publisher is required for resale or distribution outside the institutions. Permission of the publisher is required for all other derivative works, including compilations and translations. Permission of the publisher is required to store or use electronically any material contained in this journal, including any article or part of an article. For inquiries contact the publisher at the address indicated.



Lit.No: 01747 V1 02/14

Produced and published by the Scientific & Medical Affairs Department,
Smith & Nephew Inc.

Published Feb 2014

Copyright © 2014 by Smith & Nephew

Smith & Nephew Orthopaedics AG, Oberneuhofstrasse 10d, 6340 Baar,
Switzerland

◊ Trademark of Smith & Nephew. Registered US Patent and Trademark Office

Bone&JointScience is available on Smith & Nephew's Education and
Evidence website, www.smith-nephew.com/education