





# SL-PLUS<sup>◊</sup> Standard and Lateral Stem

## Table of Contents

Notes from the Author's Clinic .....	3
Indications .....	4
Contraindications.....	5
Preoperative Planning .....	5
Surgical Technique .....	6
Postoperative Treatment .....	13
Explanation of the SL-PLUS Stem.....	13
Implants.....	14
Instrumentation.....	16
Sterilization.....	20

## Nota Bene

The technique description herein is made available to the healthcare professional to illustrate the authors' suggested treatment for the uncomplicated procedure. In the final analysis, the preferred treatment is that which addresses the needs of the patient.



# Notes from the Author's Clinic

Our non-cemented stem system, which comprises a double-cone straight stem with a rectangular cross section, was implanted for the first time in 1979. Continuous follow-up resulted in ongoing improvement, culminating in the development of the SL stem in 1993. This system was suitable for the vast majority of surgical indications and provided results that were much better than those obtained with the first implants.

The SL-PLUS<sup>°</sup> Lateral shaft was added to the SL range in 2002 in order to offer additional anatomical variants which are more difficult to treat with the standard shafts in terms of optimal reconstruction of the offset.

The latest SL-PLUS straight stem that is now available integrates the proven design features of the stem large-area anchorage in the bone. However, changes have been made in the proximal section to improve intraoperative manipulation of the implant and to increase primary stability even further by improving proximal force transmission.

The extremely positive clinical results we obtained over many years, have further increased confidence in this stem design.

Professor Dr. med. K. Zweymüller  
Vienna Orthopaedic Hospital – Gersthof

## Indications Standard Stem

The SL-PLUS<sup>®</sup> Stem with and without Ti/HA (INTEGRATION) is intended for advanced hip joint wear due to degenerative, post-traumatic or rheumatoid arthritis; fracture or avascular necrosis of the femoral head.

## Indications Lateral Stem

The SL-PLUS Lateralized Stem with and without Ti/HA (INTEGRATION) is intended for varus femur forms and trumpet shape of the proximal femur (champagne flute).

These stems are for uncemented use only. These devices are intended to aid the surgeon in relieving the patient of hip pain and restoring hip motion.


## Contraindications

- Acute or chronic infections, local or systemic
- Severe muscle, nerve, or vascular diseases that endanger the respective limb
- Lack of bone substance or defective bone quality that endangers the stable seating of the prosthesis
- Any concomitant disease that may endanger the implant function
- Revision with extensive bone defects

# Preoperative Planning

Planning the correct prosthesis size is carried out preoperatively using the X-ray template. This requires an A/P X-ray and an axial X-ray. Preoperative planning should always be carried out for orientation purposes.

X-ray templates of the SL-PLUS<sup>®</sup> Standard and Lateral stem are available in a magnification of 15%.



To determine the appropriate entrypoint for access of instruments to the medullary canal, it is recommended that the surgeon draw the femoral shaft axis on the AP radiograph and extend it proximally. This line indicates how far laterally it is necessary to place the box chisel to open up the canal. This entry point is easy to locate during surgery.

After the preoperative planning, the distance from the SL-PLUS stem to the greater and lesser trochanter can be measured for the purpose of intraoperative monitoring.

For the implantation of an SL-PLUS stem, both the SL/SLR-PLUS<sup>®</sup> Basic Set as well as the corresponding SL-PLUS trial broach instrumentation are required.



# Surgical Technique

## Note

Approach: anterolateral with patient lying in the supine position.

Users of other (minimally invasive) approaches are requested to also consult the following operative instructions: posterolateral approach (Lit. No. 1426); anterior approach (Lit. No. 1494).

## Skin incision




- Previous clinical experience has shown that the SL-PLUS<sup>®</sup> stem may be successfully implanted via several skin incisions, including:
- A longitudinal incision along the anterior edge of the greater trochanter, extending from O proximally to N distally to the tip of the trochanter (Line 1).
- An oblique incision extending from the anterior edge of the greater trochanter towards the anterior superior iliac spine (Line 2).
- A reverse oblique incision approximating the intertrochanteric line (Line 3).

The fascial incision extends from the upper edge of the tip of the trochanter towards the anterior superior iliac spine.

Dorsal incision of the iliotibial band is optionally possible.

## Capsular incision and dissection




After blunt lateral entry between the tensor fasciae latae and gluteus medius/minimus, the raspatory is used to dissect along the femoral neck.

Sharp lateral retractors and a blunt medial Hohmann retractor are used during the surgical exposure. The arc of the rectus tendon is visualized, underpinned, incised, and released from its capsule.

The femoral neck is exposed via an H-shaped incision of the joint capsule, consisting of:

- A longitudinal incision, placed as far medially as possible, and extending from the acetabular margin to the intertrochanteric line and
- a proximal transverse incision of the acetabular labrum, extending around the acetabular margin from approximately the nine o'clock to the three o'clock position and
- a distal transverse incision distal extending along the intertrochanteric line.




After the wing-like opening of the joint capsule, additional dissection of the capsule can be performed by extending the distal incision in the direction of the lesser trochanter and the proximal incision medially and/or laterally.

Two blunt Hohmann retractors are positioned intra-articularly. Problematic osteophytes on the acetabular rim are removed.

The technique used for the neck resection depends on the patient (coxa vara/valga) and is selected by the surgeon (single-incision or double-incision technique).

### Capsule release




In order to facilitate alignment of the cup, an additional release of the posterior capsule is performed with the leg in the «Figure 4» position.

The knee of the operative leg is flexed, allowing it to be placed under the extended contralateral leg. In this position, the operative hip is placed in approx. 30°– 40° of adduction and 90° external rotation.

The proximal femur is mobilized with two hooks: one placed lateral to the trochanter and the other placed on the medial side of the neck of the femur.

Capsule releases have to be carried out in the direction of the lesser trochanter and the trochanteric fossa to the trochanteric tip. Additional capsular release can also be carried out on the caudal rim of the acetabulum.


## Shaft preparation



The leg is hyperextended and adducted approx. 30°–40°, and placed in 90° of external rotation, with the operative leg positioned under the opposite leg, in the «Figure 4» position.

In very muscular or obese patients, patients with a valgus femoral neck, or in cases where the proximal femur sits deep to the skin surface, further release of the posterior capsule or release of the piriformis tendon may be necessary to allow adequate mobilization of the femur prior to preparation of the implantation site.

## Entry into the medullary cavity



With the lower thigh kept in a horizontal position, the box chisel is placed close to the posterior cortex at the resection level. The box chisel has to be introduced along the femoral axis and a small square block of bone is removed. If the box chisel is not used to clear hard bone from the osteotomy site, fracture of the trochanter may occur during rasp insertion.

Driving the box chisel below the level of the resected bony surface should also be avoided.



The curved rasp facilitates opening of the diaphyseal IM canal.

Further opening of the diaphyseal medullary cavity and probing of the diaphysis with corresponding awl is recommended.



To prepare the bone bed up to size 4, start with the rasp size 01. From size 5 onwards, preparation can be started with rasp size 1. Connect the rasp with the offset adapter, which is available in a variety of different designs.

### Important


The first rasp determines the position of all the following rasps. This is why its orientation is important for an exact positioning of the stem.



Rasping is carried out using the slap hammer or the rasping machine. The weight of these instruments helps to ensure the longitudinal alignment of the rasp within the femur. It is important that a lateral force is continuously applied to the rasping machine to ensure that the rasp moves in line with the axis of the canal and does not seat in a varus position.



Attention should be paid to the anteversion and varus/valgus alignment of the rasping machine with respect to the femoral axis. Insertion of the rasps or the stem in a varus inclination increases the risk of perforation and fracture of the lateral cortex of the femur.



At the start of the rasping process, the rasp must not be inserted below the level of the estimated final position of the implant. It is extremely important to understand that the femoral osteotomy has no relationship to the final position. There is a tendency for surgeons to implant the starting smaller size rasps too deeply into the femur. This will result in an excessive enlargement of the implantation site and lead to gaps around the medial aspect of the final implant position.



The subsequent rasp is introduced into the cavity until the broach makes contact with the bone. The rasp is then driven into the femur with lateral and distal pressure using the slap hammer or the rasping machine. Then go to the next larger rasp in sequence until the appropriate size is reached.



The shoulder of the rasp corresponds to the height of the implant measured at the shoulder of the prosthesis and should correspond to the preoperatively determined distance to the greater trochanter (marked X).



The aim is to achieve a long cortical batch with a large surface area. As the IM canal is expanded in increments the contact of the cortex increases and thus also the resistance when rasping. The rasp starts to broach the edges of the shaft, i.e. the rectangular shape from the cortex. As soon as the rasp is fully working on the cortex the pitch is higher and there are signs of the yellow cortex on the edge of the rasp. The ultimate check to ensure that the size is correct is preoperative planning with the X-ray image.

In rare situations, the prosthesis size determined intraoperatively is in disagreement with the size derived from preoperative templating. If this difference is two sizes or more, the rasp may not have reached the necessary depth because of incorrect angulation or the presence of an obstacle within the canal. In such cases, the implanted prosthesis is too small to provide stable long-term fixation. In these situations, intraoperative fluoroscopy or an intraoperative radiograph should be obtained to evaluate the obstruction.



The offset adapter is removed from the detachable rasp.

### **Trial reposition**



The neck module is manually set onto the broach.



The trial ball-head can be attached to the modular neck in advance or in situ.

In each case, there is a standard modular neck for the detachable rasps of sizes 01–0, 1–6 and 7–12. The «lateral» modular necks are available to suit the detachable rasps of sizes 1–6 and 7–12.

Care should be taken that the modular neck is correctly seated on matching surface of the detachable rasp and engages properly.

The joint is repositioned and leg length, soft-tissue tension, and range of motion are checked.

If necessary the trial ball-head and/or the modular neck (standard or lateral) are changed until the results are satisfactory.

The modular neck can either be removed from the rasp manually or with a bone clamp.

The offset adapter is connected to the detachable rasp. The detachable rasp is removed from the canal using the slap hammer or the rasping machine.

## Implantation of the stem




The correct size SL-PLUS<sup>®</sup> stem is introduced manually as deep as possible into the canal, and is then seated with the impactor, using appropriately measured strokes to minimize the risk of fracture of the femur.

### Important


Pressing the stem in solely by hand is inadequate.

During impact, the protective cover remains positioned on the cone.


Once the stem is firmly seated, attempts to drive it further down the canal or to adjust its alignment within the femur cannot be performed without fracturing the bone.



Before repositioning the original ball-head the tapered trunnion cone is carefully cleaned by hand.



The ball-head is then attached to the trunnion with a slight turning motion and permanently fixed with a blow delivered with the plastic hammer.



Metal objects must never be used to deliver is used, a plastic coupling or the head impactor must be used to protect the head from direct contact with the hammer.

**Important**


Pressing the prosthetic head onto the trunnion solely by hand provides inadequate fixation.

Each femoral stem has a standard 12/14 trunnion to allow coupling with OXINIUM<sup>®</sup>, ceramic or metal heads supplied by Smith&Nephew Orthopaedics AG resp. Smith & Nephew Inc.

## Postoperative Treatment

Postoperative rehabilitation should be completed in accordance with each hospital's own practices.

## Explantation of the SL-PLUS Stem



The SL-PLUS stem can be removed using the extraction screw M8.




If this causes difficulties, the extraction block can be used.

**Important**

It is important to insert the extraction screw axially.


# Standard Implants

Item No.	SAP No.	Size	Length (L)
11420S	75002717	01	128 mm
11421S	75002719	0	132 mm
11401S	75002695	1	136 mm
11402S	75002697	2	140 mm
11403S	75002699	3	145 mm
11404S	75002701	4	150 mm
11405S	75002703	5	154 mm
11406S	75002705	6	158 mm
11407S	75002707	7	163 mm
11408S	75002709	8	168 mm
11409S	75002711	9	173 mm
11410S	75002713	10	178 mm
11411S	75002714	11	183 mm
11412S	75002715	12	188 mm




All sizes also available with Integration-PLUS<sup>®</sup> Ti-Plasma hydroxyapatite coating:

Item No.	SAP No.	Size	Length (L)
11459	75002720	01	128 mm
11460	75002722	0	132 mm
11461	75002723	1	136 mm
11462	75002724	2	140 mm
11463	75002725	3	145 mm
11464	75002726	4	150 mm
11465	75002727	5	154 mm
11466	75002728	6	158 mm
11467	75002729	7	163 mm
11468	75002730	8	168 mm
11469	75002731	9	173 mm
11470	75002732	10	178 mm
11471	75002733	11	183 mm
11472	75002734	12	188 mm




# Lateral Implants

Item No.	SAP No.	Size	Length (L)
11531S	75002748	1	136 mm
11532S	75002750	2	140 mm
11533S	75002752	3	145 mm
11534S	75002756	4	150 mm
11535S	75002758	5	154 mm
11536S	75002760	6	158 mm
11537S	75002762	7	163 mm
11538S	75002764	8	168 mm
11539S	75002766	9	173 mm
11540S	75002768	10	178 mm
11541S	75002769	11	183 mm
11542S	75002770	12	188 mm



All sizes also available with Integration-PLUS<sup>®</sup> Ti-Plasma hydroxyapatite coating:

Item No.	SAP No.	Size	Length (L)
11569	75002773	1	136 mm
11570	75002777	2	140 mm
11571	75002778	3	145 mm
11572	75002779	4	150 mm
11573	75002780	5	154 mm
11574	75002781	6	158 mm
11575	75002782	7	163 mm
11576	75002783	8	168 mm
11577	75002784	9	173 mm
11578	75002785	10	178 mm
11579	75002786	11	183 mm
11580	75002787	12	188 mm



# Instrumentation

SL-PLUS°/SLR-PLUS° Basic Instrumentation


Set No. 0942000

Item No.	SAP No.	Description	
<b>110450</b>	<b>75002198</b>	<b>Case Basic Instruments 1/4</b>	
990019	75007661	Easy Tray Lid Plastic	
① 110901	75002319	Slap Hammer	
② 110902	75002320	Extractor Block	
③ 110903	75002321	Impactor	
④ 110904	75002322	Box Chisel	
⑤ 110911	75002325	Extraction Screw M8	
⑥ 110242	75002160	Head Impactor	
⑦ 75100839*	75100839	Trial Femoral Head	22 S/+0
75100840	75100840	Trial Femoral Head	22 M/+4
75100841	75100841	Trial Femoral Head	22 L/+8
75100842*	75100842	Trial Femoral Head	22 XL/+12
75100843*	75100843	Trial Femoral Head	28 XS/-3
75100844	75100844	Trial Femoral Head	28 S/+0
75100845	75100845	Trial Femoral Head	28 M/+4
75100846	75100846	Trial Femoral Head	28 L/+8
75100847	75100847	Trial Femoral Head	28 XL/+12
75100848	75100848	Trial Femoral Head	28 XXL/+16
75100849*	75100849	Trial Femoral Head	32 XS/-3
75100850	75100850	Trial Femoral Head	32 S/+0
75100851	75100851	Trial Femoral Head	32 M/+4
75100852	75100852	Trial Femoral Head	32 L/+8
75100853	75100853	Trial Femoral Head	32 XL/+12
75100854	75100854	Trial Femoral Head	32 XXL/+16
75100855*	75100855	Trial Femoral Head	36 XS/-3
75100856	75100856	Trial Femoral Head	36 S/+0
75100857	75100857	Trial Femoral Head	36 M/+4
75100858	75100858	Trial Femoral Head	36 L/+8
75100859	75100859	Trial Femoral Head	36 XL/+12

\* special size (optional)

Optional:

Art. No.	SAP No.	Description	
75210292	75210292	SET 40 mm Trial Femoral Heads	XS/-4 to L/+8
75210293	75210293	SET 44 mm Trial Femoral Heads	XS/-4 to L/+8
⑧ 21000138	75004495	Curved Rasp	
600620	75007254	MIS Offset Box Chisel (instead of 110904)	




OPTIONAL

SL-PLUS<sup>®</sup> Trial Broach with Adapter 25 mm  
Set No. 0943001

	<b>Item No.</b>	<b>SAP No.</b>	<b>Description</b>
	<b>600930</b>	<b>75007312</b>	<b>Instrument Case</b>
	990019	75007661	Easy Tray Lid Plastic
①	600899	75007293	Trial Broach 01
	600900	75007294	Trial Broach 0
	600901	75007295	Trial Broach 1
	600902	75007296	Trial Broach 2
	600903	75007297	Trial Broach 3
	600904	75007298	Trial Broach 4
	600905	75007299	Trial Broach 5
	600906	75007300	Trial Broach 6
	600907	75007301	Trial Broach 7
	600908	75007302	Trial Broach 8
	600909	75007303	Trial Broach 9
	600910	75007304	Trial Broach 10
	600911	75007305	Trial Broach 11
	600912	75007306	Trial Broach 12
②	21000253	75004603	Trial Neck Std 01–0
	21000254	75004604	Trial Neck Std 1–6
	21000255	75004605	Trial Neck Std 7–12
	21000256	75004606	Trial Neck Lat. 1–6
	21000257	75004607	Trial Neck Lat. 7–12
③	600920	75007307	Broach Handle 25 mm
④	600923	75007310	Broach Handle Left 17/13 mm
⑤	600924	75007311	Broach Handle Right 17/13 mm

Optional:

	<b>Item No.</b>	<b>SAP No.</b>	<b>Description</b>
	600921	75007308	Adapter 40 mm (Set No. 0943005)
	600922	75007309	Adapter 10 mm (Set No. 0943000)



# Sterilization

## **Implants**

All the implants described in this Surgical Technique are sterile when they are delivered by the manufacturer. Resterilization is not allowed.

## **Instruments**

System components and instruments are not sterile when they are delivered. Before use, they must be cleaned by the usual methods in accordance with internal hospital regulations and sterilized in an autoclave in accordance with the legal regulations and guidelines applicable in the relevant country. (For detailed information please refer to leaflet Lit. No. 1363.)

The correct settings are given in the instructions for use issued by the autoclave manufacturer. Instrument manufacturers and dealers accept no responsibility for sterilization of products by the customer.



**Manufacturer**

Smith & Nephew Orthopaedics AG  
Oberneuhofstrasse 10d  
6340 Baar  
Switzerland

Telephone: +41 (0) 41 766 22 22  
Fax: +41 (0) 41 766 22 96

**Contact**

Smith & Nephew Reconstruction & Trauma  
Smith & Nephew, Inc.  
1450 Brooks Road  
Memphis, TN 38116  
USA  
Telephone: 1-901-396-2121  
Information: 1-800-8215-700  
orders and Inquiries: 1-800-238-7538