

JOURNEY[◇] II

Total Knee System


LEGION[◇]

Total Knee System

Femoral sizing guide

A/P

JOURNEY II		LEGION	
Size	(mm)	Size	(mm)
1	51.7	2	50.0
2	53.7	3N	54.0
3	56.7	3	55.0
4	59.7	4N	58.0
5	62.7	4	59.0
6	65.7	5N	61.0
7	68.8	5	62.0
8	71.8	6N	65.0
9	75.8	6	66.0
10	79.8	7	70.0
		8	75.0




◆ JOURNEY II Knee System

◆ LEGION Knee System

M/L

JOURNEY II		LEGION	
Size	(mm)	Size	(mm)
1	59.0	2	58.0
2	60.0	3N	58.0
3	61.5	3	62.0
4	64.5	4N	62.0
5	67.5	4	66.0
6	70.5	5N	66.0
7	73.5	5	70.0
8	76.0	6N	70.0
9	80.0	6	73.0
10	82.0	7	77.0
		8	80.0



◆ JOURNEY II Knee System

◆ LEGION Knee System